Surgical Procedure Introduction:

In order to help support the safety of the expanded indication for use of Vertaplex HV in the sacrum and in sacroplasty procedures and per the recommendation of the FDA, Stryker completed a cadaveric study. This study was conducted to provide data and guidelines as to the most conservative approach for sacroplasty procedures. Furthermore, this study was completed in order to provide specific guidance on performing the procedure while mitigating the potential for extravasation of cement in this new anatomical area.

Additional risks specific to the new expanded anatomical area of the sacrum are mitigated with the use of the guidance established by Stryker titled “Surgical Procedure – Long and Short Axis Sacroplasty Techniques with Post-Procedure Positioning and CT Evaluation.” The full text surgical guide will be housed for reference at www.StrykerIVS.com. A reference to this website and surgical guide will be included in the Instructions for Use.

Surgical Procedure – Long and Short Axis Sacroplasty Techniques with Post-Procedure Positioning and CT Evaluation

Technique for Long Axis Approach:

After the patient is properly placed in the prone position, the following Long axis technique is followed under (single plane) / biplane fluoroscopy (preferred):

1. Careful pre-procedural evaluation of the Sacral imaging available (MR/CT/Plain Film) that demonstrates the fracture/s and safest procedural approach, given the patients anatomical/clinical considerations, is highly recommended. A Dyna CT of the Sacrum (If available is obtained with rotational capabilities of the A Plane of the Biplane Siemens Unit. From the 3D reconstruction images, the best angle for Long Axis Approach is selected.
2. Viewing from the frontal plane, the II-A (Image Intensifier – A Plane) is rotated caudal and RAO/LAO (See #3) so the trocar access site in the Sacrum, optimizes the visualization between the SI joint laterally and the sacral foramina medially.
3. To treat the left sacral ala, the (II) image is rotated RAO so that the beam is parallel to the left sacroiliac joint. Thus, LAO for treatment of the right sacral ala.
4. Forceps are used to localize the skin insertion site at the mid-point between the inferior margin of the sacroiliac joint and the lateral margin of the left S3 neural foramen.
5. The skin and periosteum are anesthetized with lidocaine.
6. A 22-gauge spinal needle (used as a guide) is inserted into the skin and directed toward the midpoint between the superior margin of the sacroiliac joint and the lateral margin of the S1 neural foramen.

7. The II-B * If available *(Image Intensifier – B Plane) is then utilized for a lateral view for the desired needle trajectory and eventual trocar depth. Careful selection of trocar trajectory is required to avoid penetrating the anterior / posterior sacral cortex during trocar placement. Intra operative CT can also be used as an additional method for visualization during the procedure.

8. The position and orientation of the guide needle are adjusted so the needle is directed toward the center of the S1 body.

9. A skin incision is made at the base of the guide needle, and the bone cannula is inserted alongside the guide needle to the level of the periosteum.

10. The cannula is advanced into the intramedullary cavity approximately 1 cm, and the position and orientation are confirmed using frontal and lateral projections (Both II A and B planes if available).

11. The cannula is advanced using live fluoroscopy or intra operative CT on both projections until the needle tip is located in the desired position of the S1 vertebral body.

12. The cannula position is confirmed in both projections and a similar technique is followed to place the cannula in the contralateral ala. For a bilateral long axis procedure).

13. A standard mixture of Stryker® VertaPlex® HV bone cement is prepared and injected (after removal of the inner portion of the cannula / stylet) in the cannula using fluoroscopic visualization in both the frontal and lateral projections.

14. Once the desired portion of the S1 segment is adequately filled, the cannula is slowly withdrawn while continuing the cement injection until 2-4 mls of cement have been inserted.

15. When the needle approaches the inferior margin of the sacroiliac joint the cement injection is stopped, the inner cannula / stylet replaced and the cannula is removed.

Technique for Short Axis Approach:

After the patient is properly placed in the prone position, the following Short axis technique is followed under (single plane) / biplane fluoroscopy (preferred).

1. Careful pre-procedural evaluation of the Sacral imaging available (MR/CT/Plain Film) that demonstrates the fracture/s and safest procedural approach, given the patients anatomical/clinical considerations, is highly recommended. A Dyna CT of the Sacrum (If available) is obtained with rotational capabilities of the A Plane of the Biplane Siemens Unit. From the 3D reconstruction images, the best angle for Short Axis Approach will be selected.

2. Viewing from the frontal plane, the II-A (Image Intensifier – A Plane) is obliqued contralaterally from the side being injected until adequate visualization of the medial and lateral aspects of the sacroiliac joint are appreciated. This is usually obtained at a contralateral obliquity of approximately 5-25°

3. Example, to treat the right sacral ala, the image is rotated approximately 5-25° to the left (LAO) so that the beam is focused on the right sacroiliac joint. Accordingly, the image is rotated approximately 5-25 to the right (RAO) for a left sided approach.
4. Ring Forceps are used to localize the skin insertion site at the mid-point between the S1 portion of the right sacroiliac joint and the lateral margin of the right S1 neural foramen.
5. The skin and periosteum are anesthetized with lidocaine.
6. A 22-gauge spinal needle is inserted into the skin and directed toward the mid-point between the S1 portion of the right sacroiliac joint and the lateral margin of the right S1 neural foramen.
7. The II-B (Image Intensifier – B Plane) is then utilized for a lateral view for the desired needle trajectory and depth. Intra operative CT can also be used as an additional method for visualization during the procedure.
8. The position and orientation of the guide needle are adjusted so that the needle is directed toward the center of the S1 body.
9. A skin incision is made at the base of the guide needle, and the bone cannula is inserted alongside the guide needle to the level of the periosteum.
10. The cannula is advanced into S1 approximately 1 cm, and the position and orientation are confirmed using frontal and lateral projections (Both II A and B plane if available).
11. The cannula is advanced using live fluoroscopy on the lateral projection until the needle tip is located just proximal to the anterior 1/3 of the S1 segment to avoid complications within the presacral space.
12. The needle position is confirmed in both projections.
13. The same technique is repeated for advancement of a second cannula into the S2 segment, specifically the midpoint between the S2 portion of the right sacroiliac joint and the lateral margin of the right S2 neural foramen.
14. A standard mixture of Stryker® VertaPlex® HV bone cement is prepared. After removal of the inner cannula / stylet, the cement is injected in the cannulas using fluoroscopic visualization in both the frontal and lateral projections.
15. Cement is then slowly instilled via the S1 cannula under biplane live fluoroscopic guidance. Careful observation for unwanted spread of cement towards the neural foramen or through the periosteum are then performed. Cement injection can then be delayed or terminated once a desired injection pattern in the sacrum is achieved.
16. When the desired injection pattern and cement amount (2-4mls) are within the S1 level, the cement injection is stopped, the inner cannula / stylet replaced and the cannula is removed.
17. The same technique is repeated for injection and removal of the S2 cannula.

Post-Procedure Positioning:

1. Post procedure (removal of all utilized cannulas) manual compression is maintained on the puncture sites for 5 minutes to insure hemostasis
2. The patient then remains on the operative table, in the prone position, for an additional 15 minutes
3. The patient is then log rolled onto a hospital bed, where they remain supine /flat for 1 hour
4. At the one-hour point, the head of the bed may be elevated to 30° with the legs remaining extended for an additional 2 hour period.
5. 3 hours post procedure, the patient may ambulate with assistance
6. For 4-6 hours post procedure, the patient should be monitored for cardiovascular, respiratory, or neurologic dysfunction.

Post-Procedure CT Scans:

1. During the 15 minute interval that the patient remains prone on the operative table, a second Dyna CT (if available) of the Sacrum is obtained at 15 minutes with the rotational capabilities of the A Plane of the Biplane Siemens Unit. If not available, a traditional CT of the Sacrum is recommended.
2. From the CT / 3D reconstruction images, the precise location of the recently instilled sacral cement can be evaluated. Measurements / quantification can then be obtained of the desired cement trajectory within the sacrum. Additionally, occurrence, number, location, and surface area of undesired cement extravasation (if it is to occur) can be quantified and measured.

Disclaimers:

VertaPlex HV Radiopaque Bone Cement is indicated for the fixation of pathological fractures of the vertebral body using vertebroplasty or kyphoplasty. It is also indicated for the fixation of pathological fractures of the sacral vertebral body or ala using sacral vertebroplasty or sacroplasty. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).

For safe and effective use of Stryker Vertaplex HV Radiopaque Bone Cement in Sacroplasty procedures, the physician should have specific training, experience and familiarity with the injection of bone cement into weakened or diminished bone in the sacrum. Please see the Instructions for Use for the list of any contraindications for use.

This technical guide is not necessarily applicable to all patients and conditions. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product and procedure prior to its use.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient including factors such as patient's medical condition and any anatomical variations. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product

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